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The thermometer has been designed for practical use. It's not meant to replace a visit to the doctor. In order to protect your safety, please spend some times to read through this instruction and operate the device with the direction before use.

Features

Forehead Temperature

- The thermometer has been designed for practical use. It's not meant to replace a visit to the doctor. Please also remember to compare the measurement result to your regular body temperature.
- Please see the "Use of the thermometer" section to learn how to measure the body temperature.

Ear Temperature

- The thermometer has been designed for practical use. It's not meant to replace a visit to the doctor. Please also remember to compare the measurement result to your regular body temperature.
- Please see the "Use of the thermometer" section to learn how to measure the body temperature.

°C/°F Switch

In "power off" mode (The LCD only shows the room temperature). Press and hold the "SCAN" button, then immediately press and release the "ON/MEM" button. After 3 seconds, icon "°C" will be switched to icon "°F". (You may use the same way to switch back to "°C")

Scan Mode

- The Scan mode shows the actual, unadjusted surface temperatures, which is different from the body temperature. It can help you to monitor if the object temperature is suitable for the baby or patient, for example the baby's milk.
- Please see the "Use of the thermometer" section to learn how to measure the object temperature.

Fever Alarm

If thermometer detects a body temperature > 37.5°C under ear mode, there will be a long beep sound followed by three short beep sound to warn the user for potential fever.

Memory Locations

- There are total 25 sets memories for ear measurements.
- When in power on, Press the "ON/MEM" button to see the temperature stored with icon.

Last Reading

When you get a new temperature reading in ear mode, the last reading will be shown on the screen (at top right corner) with the icon.

LED Back Light

LED Back Light will automatically turned on after measurement, and automatically turned off after 5 seconds.


Room Temperature

- Suitable ambient temperature is important for the baby and patient, the thermometer always help you recognize the room temperature.
- When in power off, room temperature will always be shown on the screen with icon.
- Please see the "Use of the thermometer" section to learn how to get the correct room temperature.
- How to turn off Room Temperature: When in power on, press and hold the "ON/MEM" button for 5 sec. until "OFF" will be shown on the screen to release the button.

- Temperature measurement range: Ear/Forehead mode: 34~42.2°C (93.2~108°F), Scan mode: -22~80°C (-7.6~176°F)
- Operating temperature range: 10~40°C (50~104°F)
- Storage temperature Range: It should be stored at room temperature between -20~+50°C, RH≤85%
- Transportation temperature shall be less than 70°C, RH≤95%
- Compliance with ASTM E1965-98, EN12470-5:2003 Clinical thermometers-Part 5:Performance of infra-red ear thermometers(with maximum device), IEC/EN60601-1-2(EMC), IEC/EN60601-1(Safety) standards.
- Accuracy: For ear mode: +/-0.2°C (0.4°F) during 35.5~42°C (95.9~107.6°F), +/-0.3°C (0.5°F) for other range. For forehead mode: +/-0.2°C (0.4°F) during 36~39°C (96.8~102°F), +/-0.3°C (0.5°F) for other range. For scan mode: +/-0.3°C (0.5°F) during 22~42.2°C (71.6~108°F), others +/-4% or +/-2°C (4°F) whichever is greater.
- Fever alarm & Memory under "Ear/Forehead mode"
- °C / °F Switch function
- Battery: one lithium cell battery (CR2032 *1pcs) Ear/Forehead mode converts the ear/forehead temperature to display its "oral equivalent."
- E-LED Back Light (option) LED Back Light will automatically turned on after measurement, and automatically turned off after 5 seconds.

Specifications

Dimensions: 151.7*27.1*41.8mm
Weight (include battery): 61.1g
Power supply (Battery): DC 3V Li (CR2032X1)
Material: ABS



Important Notes

Care and cleaning	Battery replacement
<p> The probe is the most delicate part of the thermometer. Use with care when cleaning the lens to avoid damage.</p> <p>a. After the measurement, please use the cotton swab with the Alcohol (70% concentration) to clean the lens (on the inside of the probe).</p> <p>b. Allow the probe to fully dry for at least 1 minute.</p> <p>c. Storage temperature range: It should be stored at room temperature between -20~+50°C, RH≤85%</p> <p>d. Keep the unit dry and away from any liquids and direct sunlight.</p> <p>e. The Probe should not be submerging into liquids.</p> <p> * Note: Please check the device if damaged once it falls. If you can't make sure of it, please send the complete device to the nearest retailer for recalibration.</p> <p> Holding the thermometer too long may cause a higher ambient temperature reading of the probe. This could make the body temperature measurement lower than usual.</p>	<p>1. Open the battery cover: First, make sure the mode of thermometer is "forehead mode". Insert a coin to twist it clockwise, the battery is located under the twist cover at the rear of the device.</p> <p>2. Hold the device and flip the battery out with a small screwdriver</p> <p>3. Insert the new battery under the metal hook on the left side and press the right side of the battery down until you hear a click</p> <p>4. Replace the battery cover.</p> <p> Keep the battery away from children. (This device is supplied with one lithium cell (CR2032x1).) The positive (+) side up and the negative (-) side pointed down.</p> <p> Disposal of the used batteries: please comply with the local regulation.</p>


Use of the thermometer

- Always make sure the probe is clean, and without damage.
- Power on: Press "ON/MEM" button. (see figure 1)
- Mode selection:
 - 1-Ear mode: Press the "ON/MEM" button, push up the raised brim on the sliding panel on the above of LCD screen (see figure 4) until hearing a click sound, the mode of thermometer is ear mode. In this mode, please press the "SCAN" button (see figure 5) to get the body temperature by ear measurement. (see figure 6).
 - 2-Forehead mode: Press the "ON/MEM" button, push back the raised brim on the sliding panel on the above of LCD screen (see figure 6), the mode of thermometer is forehead mode. In this mode, you can get the body temperature by forehead measurement. (see figure 2)
 - 3-Scan Mode: After power on, press and hold the "ON/MEM" button, and press the "SCAN" button one time for "SCAN" mode to see icon on your LCD display. In this mode, you can get the target surface temperature. (see figure 3)

4. Temperature taking:

- 4.1-Ear temperature measuring, points for attention:
 - 4.1.1-Gently pull the ear back to straighten the ear canal (Fig.4.1.1) and snugly position the probe into the ear canal, aiming towards the membrane of the eardrum to obtain an accurate reading (Fig.4.1.2).
 - 4.1.2-Press and hold the "SCAN" button until you hear a beep sound. Remove the probe from ear and read the temperature measurement on the LCD.

Remarks:

- It is recommended that you measure 3 times with the same ear. If the 3 measurements are different, select the highest temperature.
- To avoid the risk of cross contamination, please clean the probe according to "Care and cleaning" section after each use.
- The different temperature varies in healthy persons between different parts of the body can be between 0.2~1°C
- Clinical repeatability: 0.23°C (<1 year old), 0.22°C (1~5 years old), 0.21°C (>5 years old)
- ** When you get a new temperature reading in ear mode, the last reading will be shown on the screen (at top right corner) with the icon.**

4.2-Forehead temperature measuring, points for attention:

- 4.2.1-The temporal artery is connected to the heart via the carotid artery. It is designed to measure the skin surface around temporal artery, a major artery of the head (Fig.4.2.1).
- 4.2.2-Attach the thermometer to the required measurement location on the forehead (Fig. 4.2.2). (You can choose left or right temple.)
- 4.2.3-Press the "SCAN" button, and gently scan around the temple area (Fig. 4.2.3). While scanning, you will hear a beep sound, which indicates you the newest measurement is taking place. Measurement has been completed after two short beeps are heard and the forehead icon stops flashing. The time consuming for measurement might be between 5~8 sec (up to 30sec). It depends on how much time the device needs to get the correct forehead temperature.


Remarks:

- Forehead temperature is displayed in oral mode. This mode converts the forehead temperature to display its "oral-equivalent" Value.
- Before the measurement, please stay in a stable environment for 5mins and avoid exercise, bath for 30mins.
- Remember to keep the temple area clean and away from sweat, cosmetics and scar while taking temperature.

4.3-Measuring temperature under Scan mode:

- 4.3.1-When you press the "SCAN" button, you will get the real time temperature immediately. If you press and hold the "SCAN" button, the reading of measurement will be continuous updated.
- 4.3.2-Applications include temperature measurements for Water, Milk, Cloth, Skin or other object.

Remarks:


- * Note: This mode shows the actual, unadjusted surface temperatures, which is different from the body temperature.

4.4-Get the room temperature:

- 4.4.1-For room temperature, the thermometer should be placed on the table and avoid the direct sunshine or air conditioner flow on.
- 4.4.2-When in power off, room temperature will always be shown on the screen with icon.
- 4.4.3-The room temperature can be reference after 15min later.

5. After Measurement:

- 5.1 Power off: Device will automatically shut off if left idle for more than 1 minute to extend battery life. The LCD will only show the room temperature.
- 5.2 Clean the probe after each use to ensure an accurate reading and avoid cross contamination. (See the section of Care and Cleaning for details.)
- 5.3 After measurement, please push back the raised brim on the sliding panel on the above of LCD screen for protecting the probe.

Trouble Shooting

Error Message	Problem	Solution
Er	Error 5~9, the system is not functioning properly.	Unload the battery, wait for 1 minute and repower it. If the message reappears, contact the retailer for service.
Er 1	Measurement before device stabilization.	Wait until all the icons stops flashing.
Er 2	The device showing a rapid ambient temperature change.	Allow the thermometer to rest in a room for at least 30 minutes at room temperature: 10°C and 40°C (50°F~104°F).
Er 3	The ambient temperature is not within the range between 10°C and 40°C (50°F~104°F).	
H1	1. In Ear/Forehead mode: Temperature taken is higher than +42.2°C (108°F) 2. In Scan mode: Temperature taken is higher than +80°C (176°F)	Please select the target within specifications. If a malfunction still exists, please contact the nearest retailer.
Lo	1. In Ear/Forehead mode: Temperature taken is lower than +34oC (93.2°F) 2. In Scan mode: Temperature taken is lower than -22oC (-7.6°F)	
1888	Device cannot be powered on to the ready stage.	Change with a new battery.